



BATHING RANGE TILT-IN-SPACE SHOWER CHAIR - ADULT (981-S/TIS)

PART OF THE TILT-IN-SPACE FAMILY

A fantastic Tilt-In-Space chair providing comfort to the user whilst allowing the carer ease of access during toileting and bathing routines.

FEATURES & BENEFITS

- 3D-Spacer fabric moulds to body shape providing comfort and support and fast drying
- $\bullet\,$ Recline up to 45° in-situ removing need for additional transfer
- Easy to use manual tilt mechanism with snag-free solid rod operating system requiring minimal servicing
- Safety in wet environment with braked castors
- Strong and highly durable steel main frame
- Adjustable armrests for greater ease of access
- Optional padded side-infill cushions for increased security
- Clunk-click buckle fabric attachment easily removable for washing improving hygiene and infection control

STYLE	Tilt-in-Space Chair			
SEAT WIDTH	19"(48cm)			
SWL (KG)				
OPERATION				
MATERIALS	Stainless steel frame in Satin or Powder-coated finish. 3D-Spacer fabric. Vinyl seat, arm/foot rests			
WHEEL TYPE	Attendant (Braked caster)			
HEAD SUPPORT	\checkmark			
ADJUSTABLE	Seat height, Head height, Arm rest, Foot rest			
CONFIGURABLE	\checkmark			
STANDARD COLOUR	Fabric & accessories Frame			
COLOUR OPTIONS	\checkmark			
SEAT OPTIONS	Aperture, Flexi, Full, Gap			
ARMREST OPTIONS	Straight, 'Octopus'			
FOOT REST OPTIONS	1-piece, 2-piece fold-up			
ACCESSORIES	\checkmark			
WARRANTY	15 year (frame) 2 year (components)			
COMPATIBILITY				



WARRANT	γ 15 yea	ar (frame) 2 year (com	ponents)		
COMPATIB	ILITY CLOSE	oma⊤ ∎GEBERIT			
TECHNICAL SPECIFICATIONS					
WI	DTH	LENGTH			
Seat 460mm 580mm	Overall	Overall	Overall when tilted 		
HEIGHTS					
Seat	Overall	Headrest only	Armrest		
1 560mm	1275mm	250mm	130mm		
DE	ртн	CLEARANCE	WEIGHTS		
Seat		Handle	Product	SWL	
470mm		980mm	KG 20kg	KG 180kg	

Call 033 012 80808

©2021 Care & Independence Ltd All Rig

ents of this publication are intended as a guide only. Although due care and attention has been given to ensure that instructions are acc ect, in case of doubt, please refer to your prescribing clinician or contact the Care & Independence team on 33 of 280808. Care & Inde server the right to modify, add or delete details and is under no obligation to notify the user of amendment. All data is correct at time o

